



PQ001 Oil Seeds NMR Analyzer

The PQ001 NMR Analyzer was launched in 2008. After years of upgrade, PQ001 has many advantages such as small size, high precision, good repeatability, good stability and excellent cost/benefit performance. PQ001 has been widely used in the determination of oil & moisture content and Solid Fat Content, in food & agriculture areas.




Basic parameters:

1. Magnet type: permanent magnet;
2. Magnetic field intensity: $0.25 \pm 0.05T$;
3. Probe coil: $\varnothing 40mm$;



BENCHTOP NMR FOR OIL SEEDS

- Measuring the oil content in seeds and nuts are not only essential for breeders, growers and buyers to determine their quality and commercial values, but also of importance for food manufacturers to optimize their processing.
- Niumag PQ001 NMR oil seed Analyzer provide you with a simple, fast, precise and non-destructive NMR alternative to traditional chemical techniques for oil content measurement of a variety of seeds and nuts.



NIUMAG

Niumag Corporation
 +86 15618809683 info@niumag.com
 www.nmranalyzer.com
 No,78, Keling Rd Suzhou,China